


21.3" high-bright 3 megapixel color LCDs ideal for color and grayscale diagnostic imaging applications

NEC MultiSync® MD210C3/MD211C3



21.3" color 3-Megapixel Medical Diagnostic Displays with LED Backlighting

Designed exclusively for the demanding needs of medical imaging and PACS, the NEC MultiSync MD210C3 and MD211C3, a 21.3" 3-Megapixel (MP) color display, features the flexibility of a color display with image performance rivaling that of grayscale displays.



- NEC's 3MP (1536 x 2048) LCD technology offers high brightness without compromising contrast or viewing angles, making the image quality outstanding for color or grayscale images
- 14-bit RGB lookup tables (LUTs) for gamma provide for more finely detailed, high-definition rendering of color images and crisper display of even the most delicate shadings and color differences
- LED backlight reduces startup time and lowers total power consumption while the human sensor saves power by turning off backlight when the display is not in use FDA 510(k) cleared

MultiSync MD Series Features


Factory DICOM Calibration

Every MultiSync MD Series display includes a lookup table (LUT) and is individually calibrated at the production factory to the DICOM grayscale display function for luminance. This calibration is held stable by sensors in the display for accurate out of the box viewing.




GammaCompMD™ QA software

This software is included with each display and ensures consistent image quality, including conformance to the DICOM standard, while providing an easy-to-use QA environment for medical imaging. The optional GammaCompMD QA Server provides computer networks with centralized control and management of multiple display systems.




Front Sensor

The built-in front sensor constantly monitors and maintains brightness for optimal DICOM GSDF calibration and consistent image quality. The front sensor is capable of measuring monitor brightness, white-point and contrast response, allowing it to be used for scheduled conformance, calibration and reporting functions.



ColorComp™ Digital Uniformity Control

By nature, LCD panels contain uniformity errors, which are visible as slightly brighter or darker areas on the screen. To combat this inherent trait, each MD Series display is individually characterized during production and digital uniformity correction is applied. This technology, called ColorComp, applies a digital correction to each pixel on the screen to compensate for differences in color and luminance and reduces the non-uniformity to virtually unnoticeable levels.



MODEL		MD210C3	MD211C3
DISPLAY	Viewable Image Size	21.3"	
	Native Resolution	2048 x 1536	
	Pixel Pitch	0.21mm	
	Pixels Per Inch	120 @ native resolution	
	Brightness (typical)	400 cd/m² calibrated / 740 cd/m² max.	400 cd/m² calibrated / 800 cd/m² max.
	Response Time (typical)	40ms	
	Lookup Table (LUT)	14-bit	
	Displayable Colors	16.8 million colors out of a 4.398 trillion color pallette (8-bit)	1.074 billion colors out of a 4.398 trillion color palette (10-bit) / 1024 shades of gray out of 16369 levels of grayscale
	Input Connectors	DVI-D, DisplayPort	
POWER CONSUMPTION	On (typical)	85W	
	Power Savings Mode (ADVANCED)	<2W	
PHYSICAL SPECIFICATIONS	Dimensions (WxHxD)		
	Net (with stand)	14.7 x 19.3-23.4 x 9.3 in. / 373.4 x 490.6-593.4 x 235.5mm	
	Net (without stand)	18.6 x 14.7 x 4.1 in. / 473 x 373.4 x 104.1mm	
	Weight		
	Net (with stand)	26 lbs. / 11.8 kg	
	Net (without stand)	17 lbs. / 7.8 kg	
	VESA Hole Configuration	100 x 100mm	
ENVIRONMENTAL CONDITIONS	Operating Temperature	41-95°F / 5-35°C	
	Operating Humidity	20-80%	
	Operating Altitude	10,000 ft. / 3048m	
	Storage Temperature	14-140°F / -10-60°C	
	Storage Humidity	10-85%	
	Storage Altitude	40,000 ft. / 12,192m	
LIMITED WARRANTY (parts & labor, including backlight)	5 years, including Advanced Overnight Exchange		
ADDITIONAL FEATURES	DICOM GSDF calibrated; GammaComp MD QA software; Digital CableComp; Tilt; Swivel; Height-adjustable stand, Quick Screen QA		
OPTIONAL ACCESSORIES	Color calibration sensor (MDSVSENSOR3), see http://necdisplay.com for supported video cards		
SHIPS WITH	Power cord; DVI-D cable; DisplayPort; USB cables, Quick Start Guide		

MultiSync is a registered trademark, and ColorComp and GammaCompMD QA are trademarks of NEC Display Solutions. All other brand or product names are trademarks or registered trademarks of their respective holders. Product specifications subject to change.