NEC

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INSTALLATION & MAINTENANCE GUIDE

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MD21GS-3MP

Index

Warning	Japanese.1
Contents	Japanese.2
Quick Start	Japanese.3
Display Adjustment & Mounting	Japanese.4
User Controls	Japanese.5
On-Screen Display	Japanese.6
Advanced On-Screen Display	Japanese.8
Self Calibration	Japanese.10
White Point Matching/Copy Calibration	Japanese.12
Recommended Use	Japanese.14
Specifications	Japanese.15
Features	Japanese.16
Troubleshooting	Japanese.17
References	Japanese.18

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Index

Warning	Chinese.1
Contents	Chinese.2
Quick Start	Chinese.3
Display Adjustment & Mounting	Chinese.4
User Controls	Chinese.5
On-Screen Display	Chinese.6
Advanced On-Screen Display	Chinese.8
Self Calibration	Chinese.10
White Point Matching/Copy Calibration	Chinese.12
Recommended Use	Chinese.14
Specifications	Chinese.15
Features	Chinese.16
Troubleshooting	Chinese.17
References	Chinese.18

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English

Deutsch

Français

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WARNING

TO PREVENT FIRE OR SHOCK HAZARDS, DO NOT EXPOSE THIS UNIT TO RAIN OR MOISTURE. ALSO, DO NOT USE THIS UNIT'S POLARIZED PLUG WITH AN EXTENSION CORD RECEPTACLE OR OTHER OUTLETS UNLESS THE PRONGS CAN BE FULLY INSERTED. REFRAIN FROM OPENING THE CABINET AS THERE ARE HIGH VOLTAGE COMPONENTS INSIDE. REFER SERVICING TO QUALIFIED SERVICE PERSONNEL.

CAUTION:

CAUTION

TO REDUCE THE RISK OF ELECTRIC SHOCK, MAKE SURE POWER CORD IS UNPLUGGED FROM WALL SOCKET. TO FULLY DISENGAGE THE POWER TO THE UNIT, PLEASE DISCONNECT THE POWER CORD FROM THE AC OUTLET. DO NOT REMOVE COVER (OR BACK). NO USER SERVICEABLE PARTS INSIDE. REFER SERVICING TO QUALIFIED SERVICE PERSONNEL.

This symbol warns user that uninsulated voltage within the unit may have sufficient magnitude to cause electric shock. Therefore, it is dangerous to make any kind of contact with any part inside this unit.





This symbol alerts the user that important literature concerning the operation and maintenance of this unit has been included. Therefore, it should be read carefully in order to avoid any problems.

Power Cord Important Information

CAUTION: Please use the power cord provided with this display in accordance with the table below. If a power cord is not supplied with this equipment, please contact your supplier. For all other cases, please use a power cord that matches the AC voltage of the power outlet and has been approved by and complies with the safety standard of your particular country.

When you use this Display in North America, you should use a North America Hospital Grade power cord.

Plug Type	North America Hospital Grade	European Continental	Chinese	Japanese (for Japanese Market only)
Plug Shape	 (green dot and tag label)			
Country	U.S.A./Canada	EU (except U.K.)	China	Japan
Voltage	120	230	220	100

Intended Use

The MD21GS-3MP grayscale display is intended to be used for displaying and viewing of digital images for diagnosis by trained physicians. To guarantee the display performance as specified, it must only be used in conjunction with NEC recommended display cards.

European Customers

Unpacking, Installation and calibration of this display should only be done by authorized and trained personnel. Any installation done by a non-authorized person is done under his own risk and we accept no responsibility in any device malfunctioning.

Medical Imaging


The NEC MD21GS-3MP is designed for 3-megapixel (1536 x 2048) grayscale medical imaging. This display comes fully tuned with a setting for gamma correction that complies with the DICOM Part 14 Standard.

All brands and product names are trademarks or registered trademarks of their respective owners.


Contents

Your new NEC Display display box* should contain the following:

- MD21GS-3MP Display with tilt/swivel/pivot/height adjustable stand
- Power Cord
- Quick Reference Guide
- Installation and Maintenance Guide
- Video Signal Cable (DVI-D to DVI-D cable)
- Ferrite Core^{*2}
- CD-ROM

Power Cord

CD-ROM


Quick Reference Guide

Installation & Maintenance Guide

Ferrite Core

Video Signal Cable

For complete display card instructions, please refer to the display card installation guide.

For GammaComp MD Quality Control Software, please refer to GammaComp MD User's Guide.




* Remember to save your original box and packing material to transport or ship the display.

^{*2} For Power Cord (European continental and Chinese type)

1 Open the ferrite core and damp it on the power cord.

2 Close the ferrite core.

3 Attach the Ferrite Core to the end of Power cord.
Only AC-IN (connector) side of display.

Quick Start

Display Card Installation


1. Following your PC manufacturer's guidelines, open your computer to access the PCI, AGP or PCI E slot.
2. Place the display card into an available PCI, AGP or PCI E slot and secure all screws.
3. Replace the computer cover.

NOTE: For information on installing drivers, display modes in multimonitor mode, the setting the dip switches, please refer to the display card manufacturer's documentation.

Connecting the LCD to Your PC

Once the display card is successfully installed, follow these instructions to connect the display to your PC.

1. Shut down your computer and turn off the power.
 2. Connect the DVI-D signal cable to the connector of the display card in your system (**Figure A.1**).
NOTE: For dual displays, plug a DVI-D cables into Port 1 and Port 2 of the display card (**Figure A.1**).
For further information on display card installation, refer to the display card manual. Tighten all screws.
 3. Remove the cable cover by pressing the two tabs to unlock. Lift the cover to remove (**Figure B.1**).
 4. Connect the DVI-D signal cable to the connector on the back of the display (**Figure C.1**).
NOTE: Incorrect cable connections may result in irregular operation, damage display quality/components of LCD module, and/or shorten the module life.
 5. Connect one end of the power cord to the AC inlet on the back of the display and the other end to the power outlet (**Figure C.1**). Replace the cable cover.
NOTE: Refer to the Display Adjustment & Mounting section to check tilt, raise and lower display screen and screen rotation when you manage cables and power cord.
NOTE: Please refer to the Power Cord Important Information section of this manual for proper selection of AC power cord.
 6. Collect cables and power cord and place them under the cable clip (**Figure B.1**).
NOTE: Release the cables from the clip before removing the cable cover.
 7. Turn on the display with the power switch on the top (portrait) or left side (landscape) of the display (**Figure D.1**).
Turn on the computer.
NOTE: DO NOT switch on/off repeatedly.
 8. For further adjustments, please refer to the OSD section for a full description of OSD controls.
- NOTE:** If you have any problems, refer to the Troubleshooting section.
- NOTE:** External equipment intended for connection to signal input, signal output, or other connectors, must comply with the relevant IEC standard.


Display Adjustment & Mounting

Raise and Lower Display Screen


The display may be raised or lowered in either portrait or landscape mode. To raise or lower the screen, place hands on each side of the display and lift or lower to the desired height.

NOTE: Handle with care when raising or lowering the display screen and avoid pinching your fingers.


Screen Rotation

Before rotating, the screen must be raised to the highest level to avoid knocking the screen on the stand or pinching your fingers. To rotate the screen, place hands on each side of the display screen and turn clockwise from landscape to portrait or counterclockwise from portrait to landscape. To rotate the orientation of the OSD menu between landscape and portrait, refer to the "OSD" section, "OSD Rotation" function.

Tilt and Swivel

Grasp the left and right sides of the display screen with your hands and adjust the tilt and/or swivel as desired. **NOTE:** Handle with care when tilting the screen.

Alternative VESA LCD Mounting

1. Disconnect all cables.
2. Place hands on each side of the display and lift up to the highest position. Rotate the screen 90° counterclockwise to the landscape position.
3. Place display face down on a nonabrasive surface in the landscape position. Place the screen on at least a 29 mm platform so that the stand is parallel with the surface (**Figure S.1**).
4. Remove the two top screws connecting the display to the stand (**Figure S.2**). Turn the stand to 180° counter clockwise (you will hear two clicks). Remove the screws from the bottom (**Figure S.3**) and lift off the stand. The display is now ready for mounting in an alternate manner.
5. Reverse this process to reattach stand: tighten the two bottom screws, turn stand 180° counterclockwise (you will hear two clicks), and tighten two top screws.


Figure S.1

NOTE: Use only VESA-compatible alternative mounting method (100mm pitch).

NOTE: Handle with care when removing the display stand.

NOTE: When the display is in landscape, be sure the two round indents are at the bottom of the screen (**Figure S.4**).

CAUTION: Use the original screws (4 pcs) when mounting to avoid damage to the display and stand. To fulfill the safety requirements the display must be mounted to an arm which guarantees the necessary stability under consideration of the weight of the display. The display can only be used with an approved arm (e.g. GS mark).


Figure S.2


Figure S.3



Figure S.4

User Controls

The On-Screen Display (OSD) control buttons are located on the back of the display. To access the OSD menu, press any of the user control buttons (< , > , - , +).

English

- 1. LED** Indicates power on and signal condition. The LED can be turned on, off, or dimmed. Please refer to the Advanced OSD section/LED Brightness for further information.
Green – signal is present and power is on
Orange – no signal
 - 2. RESET** Resets the highlighted control menu to the factory setting.
 - 3. SELECT** Enters the OSD controls and OSD sub menu.
 - 4/5. ADJUST -/+** Moves the bar left/right to increase or decrease the adjustment.
 - 6/7. CONTROL </>** Moves the highlighted area left/right to select control menus.
Moves the highlighted area up/down to select one of the controls.
 - 8. EXIT** Exits the OSD Controls and OSD Main Menu.
 - 9. Sensor Port** For optional external sensor used for self calibration and copy calibration.
 - 10. On-Screen Key Guide** The key guide will automatically appear on the screen when any control button on the back is pressed. The key guide will change orientation along with the OSD.
- NOTE:** When **RESET** is pressed in the main and sub-menu, a warning window will appear allowing you to cancel the **RESET** function by pressing the EXIT button.

On-Screen Display

Brightness

 **BRIGHTNESS:** Adjusts the overall image and background screen brightness. When BRIGHTNESS is adjusting, the numerical value blinks.

Position Control (except Packed Pixel mode/10 bit)

 **LEFT/RIGHT:** Controls the Horizontal Image Position within the display area.

 **DOWN/UP:** Controls the Vertical Image Position within the display area.


White Point Selection

DEFAULT: preset factory white point selection (Clear Base or Blue Base depending on model).

COPIED WHITE: White point set using copy calibration.

NATIVE: Preset native white point of the display. The maximum brightness setting is available.

Tools 1

 **SHARPNESS:** This function digitally adjusts the display for crisp images. If the factory setting is not sufficient, the image is manually adjustable. This function works when the GRAYSCALE BIT is set 8bit.

 **EXPANSION MODE:** Sets the zoom method.

FULL: The image is expanded to 1536 x 2048 regardless of the resolution.

ASPECT: The image is expanded without changing the aspect ratio.

OFF: The image is not expanded.


CUSTOM (Resolution of 1536 x 2048 only): Zoom with center fixed. Select one of nine expansion rates. In this mode the resolution may be low and show blank areas. This mode is for use with special display cards.

 **GRAYSCALE BIT:** Select the GRAYSCALE BIT (8 or 10bit) of the display card.

8 bit – 256 shades of gray

10 bit – 1024 shades of gray

Tools 2

 **LANGUAGE:** OSD control menus are available in seven languages.

 **OSD LEFT/RIGHT:** You can choose where you would like the OSD control image to appear on your top half of screen. Selecting OSD Location allows you to manually adjust the position of the OSD control menu left or right.

 **OSD DOWN/UP:** You can choose where you would like the OSD control image to appear on your top half of screen. Selecting OSD Location allows you to manually adjust the position of the OSD control menu down or up.

 **OSD TURN OFF:** The OSD control menu will remain on screen as long as it is use. In the OSD Turn Off submenu, you can select how long the display waits after the last touch of a button to shut off the OSD control menu. The preset choices are 10-120 seconds in 5 second intervals.

 **OSD LOCK OUT:** This control locks out access to all or some OSD control functions. When attempting to activate OSD controls while in the Lock Out mode, a screen will appear indicating the OSD controls are locked out. There are three types of OSD LOCK OUT:

MODE	Key strokes to enter	BRIGHTNESS	LOCK OUT Warning
OSD LOCK OUT with no control	SELECT, >	No	Yes
OSD LOCK OUT with BRIGHTNESS	SELECT, - , <	Yes	Yes
For CUSTOM LOCK OUT option, please refer to the Advanced OSD section.			

OSD LOCK OUT with no control: To activate the OSD Lock Out function, press SELECT, then ">" key and hold down simultaneously. To deactivate the OSD Lock Out, press SELECT, then ">" key and hold down simultaneously while in the OSD menu. No controls can be adjusted while in the lock out mode.

OSD LOCK OUT with BRIGHTNESS control: To activate the OSD Lock Out function, press SELECT, then "-" and "<" key and hold down simultaneously. To deactivate the OSD Lock Out, press SELECT, then "-" and "<" key and hold down simultaneously while in the OSD menu. BRIGHTNESS can be adjusted while in the lock out mode.

 **OSD ROTATION:** To rotate OSD between landscape and portrait modes. Use the + or - key to change the orientation of the OSD and On-screen Key Guide.

On-Screen Display –continued

 **RESOLUTION NOTIFIER:** This optimal resolution is 1536 x 2048. If ON is selected, a message will appear on the screen after 30 seconds, notifying you that the resolution is not at 1536 x 2048.

 **HOT KEY:** You can adjust the brightness directly. When this function is set to ON, you can adjust the brightness with <, >, +, or –, while the OSD menu is off. The standard OSD can be accessed with the EXIT button.

 **FACTORY PRESET:** Selecting Factory Preset allows you to reset all OSD control settings (BRIGHTNESS, POSITION CONTROL, WHITE POINT SELECTION, SHARPNESS, EXPANSION MODE, OSD LEFT/RIGHT, OSD DOWN/UP, OSD TURN OFF, DISPLAY MODE) back to the factory settings. Individual settings can be reset by highlighting the control to be reset and pressing the RESET button.

Information

 **DISPLAY MODE:** Provides information about the current resolution display and technical data including the preset timing being used and the horizontal and vertical frequencies.

 **MONITOR INFO.:** Indicates the model and serial numbers of your display.

OSD Warnings

To clear the on-screen warning, press the EXIT button.

NO SIGNAL: This function gives a warning when there is no Horizontal or Vertical Sync. After power is turned on or when there is a change of input signal, the **No Signal** window will appear.

RESOLUTION NOTIFIER: This function gives a warning of use with optimized resolution. After power is turned on or when there is a change of input signal or the video signal doesn't have proper resolution, the **Resolution Notifier** window will open. This function can be disabled in the TOOLS 2 menu.

OUT OF RANGE: This function gives a recommendation of the optimized resolution and refresh rate. After the power is turned on or there is a change of input signal or the video signal doesn't have proper timing, the **Out Of Range** menu will appear.

OSD LOCK OUT: This function indicates that the OSD function has been locked by the user. This can be disabled using the CUSTOM LOCK OUT menu.

LAMP LIFETIME: At 1000/500/0 hours, the OSD warning will be displayed. The display cannot hold user set luminance levels, however, white point will be maintained. Please contact your supplier when this OSD warning is displayed.

BRIGHTNESS WARNING: This function indicates when the BRIGHTNESS setting cannot be complete.

Advanced On-Screen Display

The Advanced OSD offers more control over the standard OSD controls. When OSD is off, push the "RESET" and "EXIT" button at the same time to activate the Advanced OSD. NOTE: This menu is larger than the normal OSD. To make an adjustment, ensure that the tag is highlighted, then press "SELECT". To move to another tag, press "EXIT" and then press "<" or ">" to highlight another tag.

Tag1

Brightness Adjusts the overall image and screen background brightness using the numeric value. Press "+" or "-" to adjust. When BRIGHTNESS is adjusting, a numerical value blinks.

Tag2 *This tag is not adjustable.*

Tag3

H. Position Controls Horizontal Image Position within the display. Press "+" or "-" to adjust.

V. Position Controls Vertical Image Position within the display. Press "+" or "-" to adjust.

Tag4

Gamma Selection Allows you to manually select the brightness level of grayscale. There are five selections: NO CORRECTION, DICOM, LOG LINEAR, CUSTOM, and PROGRAMMABLE.

1. NO CORRECTION

2. DICOM: DICOM GSDF (Grayscale Standard Display Function) is possible. DICOM LUT is calculated based on the calibration data at factory for each individual display.

3. LOG LINEAR: Luminance changes log rhythmically in response to a linear change of gradation.

4. CUSTOM: The gamma value is selected from a range of 0.5 to 4.0 by 0.1 intervals.

5. PROGRAMMABLE: The brightness of grayscale can be changed by using the GammaComp MD Software.

Gamma Offset Adjusts the black level of the gamma setting. **NOTE:** DICOM is not adjustable.

Tag5

White Point Selection **DEFAULT:** preset factory white point selection (Clear Base or Blue Base depending on model).

COPIED WHITE: White point set using copy calibration.

NATIVE: Preset native white point of the display. The maximum brightness setting is available.

Tag6

Self Test Mode Displays an all white screen. Used for white point matching/copy calibration.

Lamp Lifetime Estimates the life of the lamp. Lifetime depends upon luminance calibration setting.

Temperature Displays internal display temperature and internal backlight sensor temperature.

BRIGHTNESS Displays target and current luminance levels in cd/m².

INFO

Tag7

Sharpness This function is digitally adjusts the display for crisp images. If the factory setting is not sufficient, the image is manually adjustable. This function works when the GRAYSCALE BIT is 8bit.

Expansion Mode Sets the zoom method.

FULL: The image is expanded to 1536 x 2048, regardless of the resolution.

ASPECT: The image is expanded without changing the aspect ratio.

OFF: The image is not expanded.

CUSTOM (1536x2048 only): Zoom with center fixed. Select one of nine expansion rates. In this mode the resolution may be low and show blank areas. This mode is for use with special display cards.

Grayscale Bit Select the GRAYSCALE BIT, 8 or 10bit.

Language OSD control menus are available in seven languages. Press "+" or "-" to select.

OSD Left/Right You can choose where you would like the OSD control image to appear on your top half of screen. Selecting OSD location allows you to manually adjust the position of the OSD control menu left or right.

OSD Down/UP You can choose where you would like the OSD control image to appear on your top half of screen. Selecting OSD location allows you to manually adjust the position of the OSD control menu down or up.

OSD Turn Off The OSD control menu will stay on as long as it is use. In the OSD Turn Off submenu, you can select how long the display waits after the last touch of a button to shut off the OSD control menu. The preset choices are 10-120 seconds in 5 second intervals.

Advanced On-Screen Display –continued

OSD Lock Out This control completely locks out access to all OSD control functions. When attempting to activate OSD controls while in the Lock Out mode, a screen will appear indicating the OSD controls are locked out. There are three types of OSD LOCK OUT:

MODE	Key strokes to enter	BRIGHTNESS	LOCK OUT Warning
OSD LOCK OUT with no control	SELECT, >	No	Yes
OSD LOCK OUT with BRIGHTNESS	SELECT, – , <	Yes	Yes
CUSTOM LOCK OUT	<p>7 step process to Activate:</p> <ol style="list-style-type: none"> 1. Press RESET and EXIT to enter Advanced OSD 2. Press < to reach Tag 7 3. Highlight the OSD LOCK OUT function 4. Press RESET and EXIT 5. ENABLE or DISABLE: BRIGHTNESS, WARNING 6. Press SELECT to enter the OSD LOCK OUT mode 7. Press EXIT to exit the OSD <p>To Deactivate:</p> <ol style="list-style-type: none"> 1. Without OSD menu, press RESET and EXIT to bring up the LOCK OUT warning 2. Press SELECT, SELECT, <, >, <, >, EXIT 	Selectable	Selectable

OSD LOCK OUT with no control: To activate the OSD Lock Out function, press SELECT, then ">" key and hold down simultaneously. To deactivate the OSD Lock Out, press SELECT, then ">" key and hold down simultaneously while in the OSD menu. No controls can be adjusted while in the lock out mode.

OSD LOCK OUT with BRIGHTNESS control: To activate the OSD Lock Out function, press SELECT, then "–" and "<" key and hold down simultaneously. To deactivate the OSD Lock Out, press SELECT, then "–" and "<" key and hold down simultaneously while in the OSD menu. BRIGHTNESS can be adjusted while in the lock out mode.

CUSTOM: Press RESET and EXIT to enter the CUSTOM Menu. Select ENABLE or DISABLE for BRIGHTNESS or WARNING (RESOLUTION NOTIFIER/OSD LOCK OUT).

OSD Rotation To switch OSD between Landscape and Portrait modes. Press "+" or "–" to select.

Resolution Notifier The optimal resolution is 1536 x 2048. If ON is selected, a message will appear on the screen after 30 seconds, notifying you that the resolution is not set to 1536 x 2048. Press "+" or "–" to select.

HOT KEY You can adjust the brightness directly. When this function is set to ON, you can adjust the brightness with <, >, +, or –, while the OSD menu is off. The standard OSD can be accessed with the EXIT button.

LED Brightness Controls the brightness of the LED on the display, "OFF", "MID" and "MAX".

Factory Preset Selecting Factory Preset allows you to reset all OSD control settings (BRIGHTNESS, H.Position, V.Position, Gamma Selection, Gamma Offset, WHITE POINT SELECTION, SHARPNESS, EXPANSION MODE, OSD LEFT/RIGHT, OSD DOWN/UP, OSD TURN OFF, SELF CALIBRATION, DISPLAY MODE) back to the factory settings. Individual settings can be reset by highlighting the control to be reset and pressing the RESET button.

Tag8 This tag is not adjustable.

Tag9

Display Mode Provides information about the current display resolution and technical data, including the preset timing currently being used and the horizontal and vertical frequencies. Increases or decreases the current resolution. Press "SELECT" to move the adjustment mode, and press "+" or "–" to adjust horizontal and vertical resolution.

Monitor Info. Indicates the model and serial numbers of your display.

Self Calibration

Use GammaComp MD Quality Control Software for history management (Constancy test). Before proper self calibration can be performed, display should warm-up for a minimum 30 minutes. If the USB sensor is plugged in before the minute warm-up, a warning will appear on-screen (**Figure S.3**). **NOTE:** Stand-alone calibration can only be performed using the GretagMacbeth Eye-One Display Sensor. Refer to the KEY MAP (**Figure A**) when changing settings during the calibration.

Figure S.1


Figure S.2


Figure A

KEY MAP	
< / > :	changes from setting to setting
- / + :	changes the setting selection (ie...SELF or COPY)
SELECT:	moves to the next step of calibration
EXIT:	moves back one step of calibration

1. Plug-in the USB sensor to the sensor port (**Figure S.1**).
NOTE: Sensor port can only be used for the sensor connection.
2. The calibration menu will open automatically and the sensor will begin to initialize. After initializing, the DARK-MEASUREMENT message is displayed. If the DARK-MEASUREMENT message is not displayed, please proceed to step 5.
3. Place the USB sensor face down on a non-reflective, opaque surface (ie...desk or a mousepad) to measure the black level (**Figure S.2**). Press "SELECT".
4. After the DARK-MEASUREMENT has succeeded (**Figure S.4**), press "SELECT" to continue. If DARK-MEASUREMENT fails (**Figure S.5**), restart the calibration procedure.
5. Use "+" or "-" to select SELF in the MODE selection (**Figure S.6**). Press "<" or ">" to move to the next selection. Use "+" or "-" to select DEFAULT or COPIED WHITE in the WHITE selection (**Figure S.6**). Press "<" or ">" to move the next selection.
6. Select the target GAMMA: NO CORRECTION, DICOM, LOG LINEAR, PROGRAMMABLE, CUSTOM, and SKIP (**Figure S.6**).
NOTE: If CUSTOM is selected, the Gamma Value can be manually changed in intervals of 0.1 at 0.5-4.0. If SKIP is selected, the GAMMA calibration is bypassed reducing calibration time.
7. The procedure will ask for the USB sensor to be placed on the center of the display panel (**Figure S.7**). Tilt the display panel approximately 5° backward and place the USB sensor in the center of the display panel (**Figure S.1**).
NOTE: Place the USB sensor flat against the LCD to avoid external light contamination. **DO NOT** press the calibrator against the display panel. Press "SELECT".
8. Set the target LUMINANCE (cd/m²) (**Figure S.8**).

Self Calibration –continued

9. Press "SELECT" to start the calibration. Stand-alone calibration may take several minutes depending on user settings.
10. After the CALIBRATION SUCCEEDED message appears (**Figure S.9**), press SELECT.
11. To end the calibration mode, press "EXIT".


Figure S.3

Figure S.4

Figure S.5

English

Figure S.6


Figure S.7


Figure S.8

Figure S.9

White Point Matching/Copy Calibration

Use GammaComp MD Quality Control Software for history management (Constancy test). Before proper copy calibration can be performed, displays should warm-up for a minimum 30 minutes. If the sensor is plugged in before the minute warm-up, a warning will appear on-screen (**Figure C.4**). **NOTE:** Stand-alone calibration can only be performed using the GretagMacbeth Eye-One Display Sensor. Refer to the KEY MAP (**Figure A**) when changing settings during the calibration.

Figure C.1



Figure C.2

Display A – SOURCE display of white point to copy.

Display B – TARGET display which performs a copy calibration.

Figure C.3


Figure A

KEY MAP	
< / > :	changes from setting to setting
- / + :	changes the setting selection (ie...SELF or COPY)
SELECT:	moves to the next step of calibration
EXIT:	moves back one step of calibration

1. Plug-in the USB sensor to the sensor port (**Figure C.1**).
NOTE: Sensor port can only be used for the sensor connection.
2. The calibration menu will open automatically and the USB sensor will begin to initialize. After initializing, the DARK-MEASUREMENT message is displayed. If the DARK-MEASUREMENT message is not displayed, please proceed to step 5.
3. Place the USB sensor face down on a non-reflective, opaque surface (ie...desk or a mousepad) to measure the black level (**Figure C.3**). Press "SELECT".
4. After the DARK-MEASUREMENT has succeeded (**Figure C.5**), press SELECT to continue. If dark-measurement fails (**Figure C.6**), restart calibration procedure at Step 1.
5. Display A must display a white screen, enter the Advanced OSD by pressing RESET and EXIT. Use "<" or ">" keys to move into Tag 6 and select SELF-TEST MODE.
6. Using "+" or "-", select COPY in the MODE selection (**Figure C.7**). Select ON or OFF in the SELF selection. If ON selected, SELF Calibration execute with COPY Calibration.
7. Select the target GAMMA: NO CORRECTION, DICOM, LOG LINEAR, PROGRAMMABLE and CUSTOM. Please select the same target GAMMA as Display A.
NOTE: If CUSTOM is selected, the Gamma Value can be manually changed in intervals of 0.1 at 0.5-4.0.
8. The procedure will ask for the USB sensor to be placed on the center of the display panel (**Figure C.8**). Tilt the display panel approximately 5° backward and place the USB sensor in the center of the display panel (**Figure C.1**).
NOTE: Place the USB sensor flat against the display to avoid external light contamination. **DO NOT** press the calibrator against the display panel. Press "SELECT".

White Point Matching/Copy Calibration –continued

9. Press SELECT on Display B to start reading the white point information of Display A.
10. After copying Display A information, the target LUMINANCE and target chromaticity coordinates will be stored and displayed in the OSD of Display B. The TARGET VALUE is not adjustable (**Figure C.9**).
11. Remove USB sensor from source Display A and place in the center of Display B (**Figure C.2**).
12. Press SELECT to start the copy calibration. Copy calibration may take several minutes depending on user settings.
13. Once the calibration is successfully completed, the white points of Display A and Display B should match. After the CALIBRATION SUCCEEDED message appears (**Figure C.10**), press SELECT to continue into the FINE TUNING MODE.
14. At the FINE TUNING MODE message (**Figure C.11**), press EXIT if the calibration is satisfactory. If the calibration is unsatisfactory, there are two options:
 - (A) Manually fine tune the white point using the "+" and "-" keys.
 - (B) Restart the automatic copy calibration, press SELECT and return to Step 3.
15. To end the calibration mode, press "EXIT".

Figure C.4

Figure C.5

Figure C.6


Figure C.7


Figure C.8


Figure C.9


Figure C.10



Figure C.11

Recommended Use

Safety Precautions and Maintenance

FOR OPTIMUM PERFORMANCE, PLEASE NOTE THE FOLLOWING WHEN
SETTING UP AND USING THE NEC DISPLAY:

- **DO NOT OPEN THE Display.** There are no user serviceable parts inside and opening or removing covers may expose you to dangerous shock hazards or other risks. Refer all servicing to qualified service personnel.
- Do not spill any liquids into the cabinet or use your display near water.
- Do not insert objects of any kind into the cabinet slots, as they may touch dangerous voltage points, which can be harmful or fatal or may cause electric shock, fire or equipment failure.
- Do not place any heavy objects on the power cord. Damage to the cord may cause shock or fire.
- Do not place this product on an unbalanced or unstable cart, stand or table because the display may fall, causing serious damage.
- Do not place any objects onto the display and do not use the display outdoors.
- The inside of the fluorescent tube located within the display contains mercury. Please follow the bylaws or rules of your municipality to dispose of the tube properly.
- Do not bend power cord.
- Do not use display in high temperature, humid, dusty, or oily areas.
- Do not cover vent on display.

Immediately unplug your display from the wall outlet and refer servicing to qualified service personnel under the following conditions:

- When the power supply cord or plug is damaged.
- If liquid has been spilled, or objects have fallen into the display.
- If the display has been exposed to rain or water.
- If the display has been dropped or the cabinet damaged.
- If the display does not operate normally by following operating instructions.
- If glass is broken, handle with care.
- If display or glass is broken, do not come in contact with the liquid crystal and handle with care.
- Allow adequate ventilation around the display so that heat can properly dissipate. Do not block ventilated openings or place the display near a radiator or other heat sources. Do not put anything on top of display.
- The power cable connector is the primary means of detaching the system from the power supply. The display should be installed close to a power outlet which is easily accessible.
- Handle with care when transporting. Save packaging for transporting.

CAUTION

Image Persistence

Please be aware that LCD Technology may experience a phenomenon known as Image Persistence. Image Persistence occurs when a residual or "ghost" image of a previous image remains visible on the screen. Unlike CRT monitors, LCD monitors' image persistence is not permanent, but constant images being displayed for a long period of time should be avoided. To alleviate image persistence, turn off the monitor for as long as the previous image was displayed. For example, if an image was on the monitor for one hour and a residual image remains, the monitor should be turned off for one hour to erase the image.


NOTE: As with all personal display devices, we recommend displaying moving images and using a moving screen saver at regular intervals whenever the screen is idle or turning off the display when not in use.

**CORRECT PLACEMENT AND ADJUSTMENT OF THE DISPLAY CAN REDUCE EYE, SHOULDER
AND NECK FATIGUE. CHECK THE FOLLOWING WHEN YOU POSITION THE DISPLAY:**

- For optimum results during self and copy calibrations, allow 30 minutes for warm-up.
- Adjust the display height so that the top of the screen is at or slightly below eye level. Your eyes should look slightly downward when viewing the middle of the screen.
- Position your display no closer than 16 inches and no further away than 28 inches from your eyes. The optimal distance is 20 inches.
- Rest your eyes periodically by focusing on an object at least 20 feet away. Blink often.
- Position the display at a 90° angle to windows and other light sources to minimize glare and reflections. Adjust the display tilt so that ceiling lights do not reflect on your screen.
- Clean the display surface with a lint-free, non-abrasive cloth. Case of persistent dirt, wipe with cloth permeated by water, ethanol, isopropyl-alcohol completely. Avoid using any cleaning solution or glass cleaner (ex Acid, Alkali and Acetone).
- Adjust the display's brightness control to enhance readability.
- Position whatever you are looking at most of the time (the screen or reference material) directly in front of you to minimize turning your head while you are typing.
- Avoid displaying fixed patterns on the display for long periods of time to avoid image persistence (after-image effects).

Ergonomics

To obtain the best possible ergonomics benefits, we recommend the following:

- Use the preset Size and Position controls with standard signals
- Use non-interlaced signals with a vertical refresh rate of 60Hz

Cleaning the LCD Panel

- When the liquid crystal panel is stained with dust or dirt, please wipe with soft cloth gently.
- Please do not rub the LCD panel with hard material.
- Please do not apply pressure to the LCD surface.
- Please do not use OA cleaner it will cause deterioration or discolor on the LCD surface.

Cleaning the Cabinet

- Unplug the power supply
- Gently wipe the cabinet with a soft cloth
- To clean the cabinet, dampen the cloth with a neutral detergent and water, wipe the cabinet and follow with a dry cloth.

NOTE: Many plastics are used on the surface of the cabinet. DO NOT clean with benzene, thinner, alkaline detergent, alcoholic system detergent, glass cleaner, wax, polish cleaner, soap powder, or insecticide. Do not touch rubber or vinyl to the cabinet for a long time. These types of fluids and fabrics can cause the paint to deteriorate, crack or peel.

Specifications

Display Specifications		MD21GS-3MP Display	Notes
LCD Module	Diagonal:	54 cm/21.3 inches	Active matrix; thin film transistor (TFT) liquid crystal display (LCD); 0.2115mm dot pitch; 400 cd/m ² calibrated white luminance; 700:1 contrast ratio, typical.
	Viewable Image Size:	54 cm/21.3 inches	
	Native Resolution (Pixel Count):	1536 x 2048 (Portrait) 2048 x 1536 (Landscape)	
Input Signal	Video:	Digital Input: DVI	VideoDot Clock 162Mhz Max
Grayscale Tones	10 Bit:	1024 (10bit) shades of gray from a pallet of 3061 (11.5bit)	When used with an recommended 10 bit display card
	8 Bit:	256 (8bit) shades of gray from a pallet of 3061 (11.5 bit)	When used with an recommended 8 bit display card
Synchronization Range	Horizontal:	31.5 kHz to 96.0 kHz	Automatically
	Vertical:	30 Hz, 50 Hz to 85 Hz	Automatically
Viewing Angle	Left/Right:	±88° (CR > 10)	
	Up/Down:	±88° (CR > 10)	
Image Formation Time		35 ms (Typ.)	
Active Display Area	Portrait: Horiz.:	325 mm/12.8 inches	
	Vert.:	433 mm/17.1 inches	
	Landscape: Horiz.:	433 mm/17.1 inches	
	Vert.:	325 mm/12.8 inches	
Power Supply		AC 100-240V ~ 50/60Hz	
Current Rating		0.7 - 0.35 A	
Dimensions	Portrait:	358.7 mm (W) x 514 - 560 mm (H) x 200 mm (D) 14.1 inches (W) x 20.2 - 22.0 inches (H) x 7.9 inches (D)	
	Landscape:	467 mm (W) x 460 - 506 mm (H) x 200 mm (D) 18.4 inches (W) x 18.1 - 19.9 inches (H) x 7.9 inches (D)	
	Height Adjustment:	46 mm/1.8 inches	
Weight		11.7 kg/25.8 lbs	
Environmental Considerations			
	Operating Temperature:	5°C to 40°C/41°F to 104°F	
	Humidity:	30% to 80% (without condensation)	
	Feet:	0 to 10,000 Feet/0 to 3,048 m	
	Storage Temperature (in package):	-20°C to 60°C/-4°F to 140°F	
	Humidity (in package):	10% to 85% (without condensation)	
	Feet (in package):	0 to 40,000 Feet/0 to 12,192 m	

UL-Classification

According to the type of protection against electric shock: CLASS I

According to the degree of protection against electric shock: No Patient connection

According to the degree of protection against ingress of water as detailed in the current edition of IEC529: No Protection

According to the method of sterilization or disinfection recommended by the manufacturer: Not Specified

According to the degree of safety of application in the presence of a FLAMMABLE AN AESTHETIC MIXTURE WITH AIR or a WITH OXYGEN OR NITROUS OXIDE: Not suitable

According to the mode of operation: Continuous operation

Features

DVI-D: The digital-only subset of DVI ratified by the Digital Display Working Group (DDWG) for digital connections between computers and displays. As a digital-only connector, analog support is not provided off a DVI-D connector. As a DVI-based digital only connection, only a simple adapter is necessary for compatibility between DVI-D and other DVI-based digital connectors such as DFP and P&D.

P&D (Plug and Display): The VESA standard for digital flat panel display interfaces. It is more robust than DFP since it allows for other options off a signal connector (options like USB, analog video and IEEE-1394-995). The VESA committee has recognized that DFP is a subset of P&D. As a DVI-based connector (for the digital input pins), only a simple adapter is necessary for compatibility between P&D and other DVI-based digital connector such as DVI and DFP.

Pivoting Stand: Allows users to adjust the display to the orientation that best fits their application, either Landscape orientation for wide documents, or portrait orientation for the ability to preview a full page on one screen at one time.

Reduced Footprint: Provides the ideal solution for environments requiring superior image quality but with size and weight limitations. The display's small footprint and low weight allow it to be moved or transported easily from one location to another.

OSD (On-Screen-Display) Controls: advanced OSD controls allow you to quickly and easily adjust all elements of your screen image via simple to use on-screen menus and includes a complete OSD Lock Out function.

ErgoDesign Features: Enhance human ergonomics to improve the working environment, protect the health of the user and save money. Examples include OSD controls for quick and easy image adjustments, tilt base for preferred angle of vision, small footprint and compliance with MPRII guideline for lower emissions.

Anti-glare and Low-reflection Screen: an anti-glare and low-reflection screen reduces glare and ambient reflection. The BLACK level perception of the display is improved resulting in a clearer image.

Multiple Frequency Technology: Automatically adjusts display to the display card's scanning frequency, thus displaying the resolution required.

FullScan Capability: Allows you to use the entire screen area in most resolutions, significantly expanding image size.

Wide Viewing Angle Technology: Allows the user to be able to see the display from any angle (176 degrees) from any orientation — Portrait or Landscape. Provides full 176° viewing angles either up, down, left or right.

VESA Standard Mounting Interface: Allows users to connect their display to any VESA standard third party mounting arm or bracket. Allows for the display to be mounted on a wall or an arm using any third party compliant device.

X-Light Backlight System: Using an internal sensor, this allows for consistent brightness and white point levels over a long period of time, providing the basis for excellent diagnostic quality. The system also allows for easy multiple display white point matching with the use of a calibrator puck (sold separately), contributing to reduced costs and increased efficiency.

Troubleshooting

No picture

- The signal cable should be completely connected to the display card/computer.
- The display card should be completely seated in its slot.
- Check the power switch should be in the ON position.
- Computer power switch should be in the ON position.
- Check to make sure that a supported mode has been selected on the display card or system being used.
(Please consult display card or system manual to change graphics mode.)
- Check the display and your display card with respect to compatibility and recommended settings.
- Check the signal cable connector for bent or pushed-in pins.

Power Button does not respond

- Unplug the power cord of the display from the AC outlet to turn off and reset the display.
- Check the power switch on the top of the display (portrait) or side of the display (landscape).

Image Persistence

- Please be aware that LCD Technology may experience a phenomenon known as Image Persistence. Image Persistence occurs when a residual or "ghost" image of a previous image remains visible on the screen. Unlike CRT monitors, LCD monitors' image persistence is not permanent, but constant images being displayed for a long period of time should be avoided. To alleviate image persistence, turn off the monitor for as long as the previous image was displayed. For example, if an image was on the monitor for one hour and a residual image remains, the monitor should be turned off for one hour to erase the image.

NOTE: As with all personal display devices, we recommend using a moving screen saver at regular intervals whenever the screen is idle or turning off the display when not in use.

Message "OUT OF RANGE" is displayed (screen is either blank or shows rough images only)

- Image is roughly displayed (pixels are missing) and OSD warning "OUT OF RANGE" is displayed means the signal clock or resolution is too high. Choose one of the supported modes.
- OSD warning "OUT OF RANGE" is displayed on a blank screen: Signal frequency is out of range. Choose one of the supported modes.

Image is unstable, unfocused or swimming is apparent

- Signal cable should be completely attached to the computer.
- Check the display and your display card with respect to compatibility and recommended signal timings.
- If your text is garbled, change the video mode to non-interlace and use 60Hz refresh rate.

Display image is not sized properly

- Check to make sure that a supported mode has been selected on the display card or system being used.
(Please consult display card or system manual to change graphics mode.)
- Check the EXPANSION MODE in the OSD.

No Video

- If no video is present on the screen, turn the Power button off and on again.
- Make certain the computer is not in a power-saving mode (touch the keyboard or mouse).

No Tilt

- Rotate the screen to 90° clockwise until you hear one click.
- Refer to the Display Adjustment and Mounting section to insure correct installation of stand.

NO OSD

- Refer to the OSD LOCK OUT function in the Advanced OSD section of this manual.

References

North America and Canada

NEC Monitor Customer Service & Support

Customer Service and Technical Support:

(800) 632-4662

Fax: (801) 907-3805

Parts and Accessories/Macintosh

Cable Adapter:

(800) 632-4662

Warranty Information:

www.necddisplay.com

Online Technical Support

www.necddisplay.com

Sales and Product Information

Sales Information Line:

(888) 632-6487

Canadian Customers:

(866) 771-0266, Ext#: 4037

Government Sales:

(800) 284-6320

Government Sales email:

gov@necddisplay.com

Electronic Channels

World Wide Web:

www.necddisplay.com

Product Registration:

www.necddisplay.com

Drivers and Downloads

www.necddisplay.com

Europe

Address

NEC Display Solutions Europe GmbH

Phone: +49(0)89 99699 - 0

Landshuter Allee 12 -14

D-80637 München

Electronic Channels

World Wide Web:

www.medical.nec-display-solutions.com

DECLARATION OF CONFORMITY

For USA

This device complies with Part 15 of FCC Rules. Operation is subject to the following two conditions. (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

U.S. Responsible Party:	NEC Display Solutions of America, Inc.
Address:	500 Park Blvd, Suite 1100
	Itasca, Illinois 60143
Tel. No.:	(630) 467-3000

Type of Product:	Display Monitor
Equipment Classification:	Class B Peripheral
Model:	MD21GS-3MP-CB/MD21GS-3MP-BB

*We hereby declare that the equipment specified above
conforms to the technical standards as specified in the FCC Rules.*

FCC Information

- Use the attached specified cables with the MD21GS-3MP grayscale display so as not to interfere with radio and television reception.
 - Please use the supplied power cord or equivalent to ensure FCC compliance.
 - Please use the supplied shielded video signal cable.Use of other cables and adapters may cause interference with radio and television reception.
- This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy, and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:
 - Reorient or relocate the receiving antenna.
 - Increase the separation between the equipment and receiver.
 - Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
 - Consult your dealer or an experienced radio/TV technician for help.

If necessary, the user should contact the dealer or an experienced radio/television technician for additional suggestions. The user may find the following booklet, prepared by the Federal Communications Commission, helpful: "How to Identify and Resolve Radio-TV Interference Problems." This booklet is available from the U.S. Government Printing Office, Washington, D.C., 20402, Stock No. 004-000-00345-4.

For Europe

Means of Conformity

Device Classification:	Class I, non-measuring function
Applicable Rules:	Annex IX, Rules 1.4 (Section 1) and 1.1 (Section 3)
Product Name:	NEC MD21GS-3MP 21.3" Diagnostic Imaging Grayscale LCD Monitor
Model Numbers:	MD21GS-3MP-BK-BB (Blue Base) MD21GS-3MP-BK-CB (Clear Base)
UMDNS Code:	16-603

NEC Display Solutions Europe GmbH declares that the products listed are in conformity with the essential requirements and provisions of the Council Directive 93/42/EEC and conform to the applicable clauses of the following standards:

– EN 60601-1
– EN 60601-1-2
– EN 61000-3-2
– EN 61000-3-3

NEC Display Solutions Europe GmbH
Landshuter Allee 12-14. 80637 München, Germany

For Canada

Canadian Department of Communications Compliance Statement

DOC: This Class B digital apparatus meets all requirements of the Canadian Interference-Causing Equipment Regulations.
C-UL: Bears the C-UL Mark and is in compliance with Canadian Safety Regulations according to CAN/CSA C22.2 No. 601.1.

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